Glidewire GT-R™

Guide Wire / Guide

ENGLISH

Please read all instructions prior to use.

CAUTION: Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

INDICATIONS FOR USE

The Glidewire GT-R is intended to direct a catheter to the desired anatomical location in the peripheral vasculature during diagnostic or interventional procedures. This device is not intended for neurovascular or coronary interventions.

DESCRIPTION

The Glidewire GT-R is a wire made of super-elastic Nickel-Titanium alloy wire with a coiled radiopaque tip marker which is entirely coated with a polyurethane layer containing tungsten (inner layer), and secondly coated with a hydrophilic polymer (outer layer). Due to the presence of the hydrophilic coating, the surface of the guide wire becomes lubricious once wet with saline solution or blood and thereby reduces friction resistance.

Product	Min. Compatible microcatheter I.D.
Glidewire GT-R 0.016"	0.017" (0.43mm)
Glidewire GT-R 0.018"	0.022" (0.57mm)

CONTRAINDICATIONS

No-contraindications are known when using this device in peripheral vascular procedures.

Carefully read the list of Contraindications in the instructions for use accompanying the other interventional devices to be used with the guide wire.

COMPLICATIONS

Complications that may be related to the use of guide wire can include but may not be limited to the following :

• vessel perforation • intimal tear • vessel damage • vessel spasm • hematoma at the puncture site • infections

WARNINGS/PRECAUTIONS

Failure to abide by the following warnings might result in damage to the vessel, shearing of the Glidewire GT-R, and damage to the device coating (referring to hydrophilic coating and polyurethane jacket), which may necessitate intervention or result in serious adverse events.

WARNINGS

- Do not manipulate or withdraw the Glidewire GT-R through a metal entry needle, a metal dilator, or a metal guide wire inserter. Manipulation and/or withdrawal through a metal entry needle or a metal dilator may result in destruction and/or separation of the device coating requiring retrieval. A plastic entry needle is recommended when using this wire for initial placement.
- Do not use the Glidewire GT-R with a metal-tipped catheter, a laser catheter, or a catheter that has a rotating cutting design such as an atherectomy catheter. Use of the Glidewire GT-R with such catheters may result damage to the surface of the wire and/or separation of the wire.
- Do not attempt to reshape the Glidewire GT-R unless it is indicated as "shapeable" on the package. Attempting to reshape the wire may cause damage, resulting in the release of wire fragments into the vessel.
- When reshaping the shapeable Glidewire GT-R tip use mandrel. Do not use heat, sharp metal instruments, or excessive force. These methods may result in abrasion of hydrophilic coating and damage to the wire.
- When exchanging or withdrawing a catheter over the Glidewire GT-R, secure and maintain the guide wire in place under fluoroscopy to avoid unexpected guide wire advancement, otherwise damage to the vessel wall by the wire's tip may occur.
- A retrieving device, such as a gripper or basket forceps, can only be used after the Glidewire GT-R has been removed from the patient's vessel. Using a retrieving device while the Glidewire GT-R is in the vessel may cause the Glidewire GT-R to break.
- Manipulate the Glidewire GT-R slowly and carefully in the vessel while confirming the behavior and location of the wire's tip under fluoroscopy to avoid damage to the vessel. If any resistance is felt, the behavior and/or position of the wire's tip seems improper, the wire is kinked, or vessel spasm is suspected, stop manipulating the wire (and the catheter) and determine the cause carefully by fluoroscopy to take suitable remedial actions. Then remove the wire slowly without turning. Failure to exercise proper caution may result in bending, kinking, separation of the guide wire's tip, damage to the catheter, or damage to the vessel.
- Do not turn the proximal end of the Glidewire GT-R three or more turns in succession in the same direction if the wire's tip is trapped due to vessel spasm or some other cause. Separation of the wire may result.
- Carefully manipulate the guiding catheter and observe it for undesirable behaviors while the Glidewire GT-R is in the vessel. Especially when the wire's tip is trapped due to vessel spasm or some other cause, damage to or separation of the wire and/ or damage to the vessel may result if the guiding catheter slips out of position and exerts too large a load on the wire.
- Do not attempt to use the Glidewire GT-R if it has been bent, kinked or damaged. Use of a damaged wire may result in damage to the vessel or the release of wire fragments into the vessel.
- Consider the use of systemic heparinization to prevent or reduce the possibility of thrombus formation on the surface of the Glidewire GT-R.
- The Glidewire GT-R may slide entirely into the catheter or slip out of the catheter because of its low sliding friction.

PRECAUTIONS

- The Glidewire GT-R is intended to be used by a physician, who is well trained in manipulation and observation of guide wires under fluoroscopy.
- The Glidewire GT-R has been sterilized by ethylene oxide gas. For single use only. Do not reuse. Do not resterilize. Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.
- Sterile, non-toxic and non-pyrogenic in unopened and undamaged unit package. Carefully inspect the unit package to make sure that no damage to the sterile barrier has occurred. Do not use if the unit package or the guide wire is broken or soiled.
- The Glidewire GT-R should be used immediately after opening the unit package, and be disposed of safely and properly after use following the local regulations for medical waste management.
- The entire operation should be carried out aseptically.
- When using a drug or a device concurrently with the Glidewire GT-R, the operator should have a full understanding of the properties/characteristics of the drug or device so as to avoid damage to the Glidewire GT-R.
- The surface of the Glidewire GT-R is not lubricous unless it is wet. Before taking it out of its holder and inserting it through a catheter, fill the holder and the catheter with heparinized saline solution.
- When reinserting the Glidewire GT-R back into the holder, take care not to damage the wire's hydrophilic coating with the edge of the holder.
- Use the enclosed RADIFOCUS TORQUE DEVICE to handle the Glidewire GT-R. Do not use a metal torque device with the Glidewire GT-R. Use of a metal torque device may result in damage to the wire. Also, do not slip a tightened up torque device or Y-connector over the wire, as this may result in damage to the wire.
- Due to variations of certain catheter tip inner diameters, abrasion of the hydrophilic coating may occur during manipulation. If any resistance is felt during introduction of the catheter, it is advisable to stop using such catheters.
- Using an incompatible device with the Glidewire GT-R may cause damage to the hydrophilic coating, and/or damage the wire.
- Do not manipulate the Glidewire GT-R through a tightened up rotating hemostasis valve, as this may result in damage to the wire.
- After removal from the patient's vessel, and prior to reinserting it into the same patient during the same catheterization, the Glidewire GT-R should be rinsed in a bowl full of heparinized saline solution. Any blood residues still adhering to the wire can be removed by wiping once with a gauze moistened with heparinized saline solution. Avoid excessive wiping or wiping with dry gauze as this may damage the coating of the guide wire. Use of alcohol, antiseptic solutions or other solvents must be avoided, because they may adversely affect the surface of the Glidewire GT-R.
- Do not use under MRI. This device is not MRI compatible.
- This device is intended to be used under fluoroscopic guidance, determine if therapeutic benefits outweigh risks in individuals who are or may be pregnant.
- This device is to be used under fluoroscopy, all protective measures regarding radiation are to be respected.
- The safety and effectiveness of the device is unknown in vascular regions other than those specifically indicated.

DIRECTIONS FOR USE

- Remove the Glidewire GT-R and the holder together from the package.
- Fill the holder with heparinized saline solution through the Wire clip or hub of the holder using a syringe.
- Remove the Protective cap and Wire clip from the holder.
- Remove the inserter and then the Glidewire GT-R from the holder, and inspect the wire prior to use to verify it is lubricated. If the Glidewire GT-R can not be easily removed from the holder, inject more heparinized saline solution into the holder and try again.
- Carefully reshape the tip of the shapeable wire by using the mandrel while the surface is wet.
- Prior to use prime the catheter with heparinized saline solution to ensure smooth movement of the Glidewire GT-R within it.
- Keep at least 5 cm of the wire extended out of the hub of the catheter during introduction. A RADIFOCUS TORQUE DEVICE, which is included in the package, is recommended for easier handling/manipulation of the wire.
- When using the Torque Device, loosen the nut and slip it over the proximal end of the Glidewire GT-R.
- When the Torque Device is in the desired position on the Glidewire GT-R, tighten the nut. The Torque Device may be repositioned or removed by loosening the nut and retightening it.
- Manipulate the Glidewire GT-R in the vessel while monitoring the behavior and location of the wire's tip under fluoroscopy.

PRECAUTIONS FOR STORAGE

Avoid exposure to water, direct sunlight, extreme temperatures and high humidity during storage. Store under controlled room temperature.

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