**Pinnacle®**  
**Pinnacle® Precision**  
**Access System®**  
**Pinnacle® TIF Tip™**

Introducer Sheath / Gaine d'introduction  
Vaina Introductora / Bainha Introdutora

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TERUMO and Pinnacle are registered trademarks of TERUMO CORPORATION.  
PRECISION ACCESS SYSTEM is a registered trademark of TERUMO MEDICAL CORPORATION.

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## DESCRIPTION AND INDICATIONS FOR USE

The PINNACLE INTRODUCER consists of an introducer (a sheath and a dilator), a mini guidewire, a guide inserter, and may contain an introducer needle.

The Pinnacle Introducer is used to facilitate placing a catheter through the skin into a vein or artery. The mini guidewire is an accessory device which is used for placement of the sheath into the vein or artery. The introducer needle is an accessory device which is used in facilitating entry through the skin into a vein or artery to provide a conduit for the mini guidewire.

## COMPONENT DESCRIPTION

Refer to product labeling for appropriate system components. All components that enter blood vessels are radiopaque.

### Sheath

Incorporates a 1-way valve and a 3-way stopcock connected by a side tube. The sheath can be used with a catheter of the same Fr. size or up to two Fr. sizes smaller without blood leakage at the 1-way valve. This highly flexible sheath is designed to resist kinking.

### Dilator

The precise fit of the dilator in the sheath allows for simultaneous motion of both dilator and sheath.

### Spring Coil Mini Guidewire or Nitinol Mandrel Guidewire


A 45 cm stainless steel guidewire or 43 cm Nitinol guidewire with 0.21 inch, 0.35 inch or 0.38 inch outside diameter is included depending on the catalog number of the system.

### Guide Inserter

Before inserting the mini guide wire, the guide inserter can be set at the entry of the needle hub for easy insertion.

### Needle

The introducer needle is an accessory device which is used in facilitating the entry through the skin into a vein or artery to provide a conduit for the mini guidewire. The needle has a tapered distal end to provide 21Ga access to vessel. Needle may be included depending on catalog number of the system.

The contents of this package may differ from that shown in the instructions for use. Please consult available catalog information for proper contents.

## PRECAUTIONS

- When using metal needle cannula, do not withdraw the guidewire back into the cannula, as shearing of the guidewire may result.
- This kit must be used by a trained physician.
- This kit is for single use only. Do not resterilize or reuse.

- Contents are sterile, non-toxic, and non-pyrogenic in unopened, undamaged package.
- Do not use if the package or product is stained or damaged. Use the introducer kit immediately after opening the package and dispose of the kit after use.
- Before use, make sure the sheath size (Fr.) is appropriate for the access vessel and the catheter to be used.
- The entire procedure from skin incision to sheath removal must be carried out aseptically.
- Do not use a power injector through the side tube and 3-way stopcock.

**Caution**


- Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

## INSTRUCTIONS FOR USE

1. If necessary, make a small skin incision at the puncture site with a surgical knife.
2. Insert a cannula into the vessel (Fig. 1).
3. Insert the selected flexible end of the mini guidewire through the cannula into the vessel (Fig. 2).

**Caution**

- Advance the mini guidewire slowly. If resistance is met, do not advance or withdraw the mini guidewire until the cause of resistance is determined.
4. Remove the cannula over the mini guidewire.
  5. Connect a flushing line to the 3-way stopcock of the introducer sheath. Fill the sheath assembly completely with heparinized saline, removing all air.
  6. Prime the dilator using a syringe with heparinized saline.
  7. Insert the vessel dilator fully into the sheath. The female hub of the sheath connects with the male hub of the dilator, and locks in place by means of grip (Fig. 3).





**Cautions**

- Insert the dilator into the center of the sheath valve. Forced insertion of the dilator which misses the center of the sheath valve may cause damage, and result in blood leakage.
  - Lock the dilator hub into the sheath hub securely. If the dilator hub is not locked into the sheath hub, only the sheath will advance into the vessel and the tip of the sheath may damage the vessel. Advancing the sheath alone may cause damage to the vessel.
8. Insert the dilator and sheath together over the mini guidewire, and into the blood vessel (Fig. 4).
  9. "Unlock" the dilator hub from the sheath hub by bending the dilator hub downward (Fig. 5).
  10. Slowly remove the dilator and mini guidewire together, leaving the sheath in the vessel (Fig. 6). If injection or sampling is necessary at this point, remove the mini guidewire only, and use the dilator hub as an injection port before removing it.

ENGLISH

**Caution**

- Slowly remove the dilator from the sheath. Rapid withdrawal of the dilator may result in the incomplete closing of the 1-way valve, resulting in blood flow through the valve. If this occurs, replace the dilator into the sheath and remove again slowly.
11. Insert a catheter through the sheath and into the blood vessel, and advance to the desired location (Fig. 7).

### Cautions



- Before removing or inserting the catheter through the sheath, aspirate blood from the 3-way stopcock to remove any fibrin deposition which may have accumulated in or on the tip of the sheath (Fig. 8).
- When puncturing, suturing, or incising the tissue near the sheath, be careful not to damage the sheath. Do not put a clamp on the sheath nor bind it with a thread.

The flushing line may also be used as a continuous infusion site by connecting an infusion line to the 3-way stopcock (Fig. 9).

12. When inserting, manipulating or withdrawing a catheter from the sheath, always hold it in place. To temporarily suture the sheath (for continued access) use the suture eye.

### Cautions

- Do not place suture on the sheath tubing since this may restrict access/flow through the sheath.
  - Do not use a power injector through the side tube and 3-way stopcock.
13. When exchanging catheters, remove the used catheter and repeat Step 11.
  14. After the intended procedure is completed, remove the catheter and then the sheath.

### CAUTION FOR STORAGE

Do not store at extreme temperatures and humidity. Avoid direct sunlight.

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