




***Glidewire GT-R™***

Guidewire

# RESHAPING PRECISION

**Designed with enhanced shapeability<sup>1-2</sup>**

- Soft, flexible tip to track through visceral anatomy with ease
- Superior tip shape retention
- Kink-resistant, nickel-titanium core
- Available in up to 220cm for radial access

---

***0.016" and 0.018" diameters compatible with  
the PROGREAT® family of microcatheters***

---

 **TERUMO**  
INTERVENTIONAL  
SYSTEMS

# PUSHING BOUNDARIES

Terumo Interventional Systems is **committed to your success** with innovative procedural solutions and ongoing support for your most challenging cases.

We are relentlessly seeking new ways to help you achieve **better outcomes for more patients.**

## GLIDEWIRE GT-R™ Guidewire

| PRODUCT CODE | DIAMETER | LENGTH | TIP SHAPE              | TAPER LENGTH | RADIOPAQUE LENGTH |
|--------------|----------|--------|------------------------|--------------|-------------------|
| GT*ES1618A   | 0.016"   | 180 cm | Straight/<br>Shapeable | 25 cm        | 3.5 cm            |
| GT*ES1818A   | 0.018"   | 180 cm |                        | 30 cm        |                   |
| GT*ES1622A   | 0.016"   | 220 cm |                        | 25 cm        |                   |
| GT*ES1822A   | 0.018"   | 220 cm |                        | 30 cm        |                   |

Superior hydrophillic coating



3.5cm gold coil  
at distal tip  
for optimal visibility

## FIND OUT MORE



US: 800.888.3786  
CA: 833.883.7866



US: [terumo.com](http://terumo.com)  
CA: [terumocanada.ca](http://terumocanada.ca)

**RX ONLY.** Refer to the product labels and package insert for complete warnings, precautions, potential complications, and instructions for use.

### References:

1. Data on file.
2. Glidewire GT-R 510(k).

©2025 Terumo Medical Corporation and Terumo Medical Canada Inc. All rights reserved.  
All brand names are trademarks or registered trademarks of their respective owners. PM-09548

 **TERUMO**  
INTERVENTIONAL  
SYSTEMS